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
Detection of- histidine tagged fusion proteins by using a- high specific mouse monoclonal- anti histidine tag antibody

© by Zentgraf, H.1, Frey, M.1, Schwinn, S.1, Tessmer, C.1, Willemann, B.1, Samstag, Y.21, Velhagen, I.1 Nucl. Acids Res., 23, 3347-3348 (1995).

Recombinant proteins are utilized for various purposes in molecular biology such as the production of antibodies or investigation of the mechanism of protein-nucleic acid or protein-protein interactions. Many prokaryotic expression vectors have been established enabling synthesis of the protein of interest as a fusion with a peptide, thus facilitating purification. Expression of recombinant proteins in *Escherichia coli* as a fusion protein with neighbouring histidine residues is one of the most popular methods (1; Qiagen), because these proteins have useful attributes. The affinity of the histidine tag motif to Ni²⁺ by chelation is strong and selective enough to enable purification of the protein to homogeneity by affinity chromatography on a Ni²⁺-NTA adsorbant (1). Additionally, an enterokinase cleavage site engineered between the N- or C- terminus of the protein of interest and the polyhistidine fusion partner facilitates the removal of the tag. A drawback of the polyhistidine fusion system so far may be the lack of sera specific for the affinity tag. To circumvent this limitation we raised mouse monoclonal anti-histidine tag antibodies by immunizing mice with mixtures histidine-tagged proteins, engineered in our laboratory using standard immunizing protocols (2). Such an antibody allows the detection of proteins with a histidine tag without the need to raise specific sera against the protein of interest. Here we present data obtained using the anti-histidine tag antibody mAb 13/45/31(subclass IgG2a).


Fig. 1 Immunofluorescence of SAOS-2 cells transfected with the expression plasmid pUHD-His-p53mt267. Cells were stained simultaneously with (a) DAPI, (b) mAb 13/45/31 as first and a FITC coupled goat-anti-mouse antibody as secondary antibody and (c) a polyclonal rabbit anti-p53 serum and, subsequently, a Cy3 coupled goat-anti-rabbit antibody.

First we tested mAb 13/45/31 in immunofluorescence analysis. For this we transfected p53-negative SAOS-2 cells with the expression plasmid pUHD-His-p53mt267. This plasmid is a derivative of plasmid pUHD 10-1 (3) with additional six histidine residues and an enterokinase cleavage site at the N-terminus of the p53 coding region. Staining with DAPI (4',6-diamidine-2'-phenylindole dihydrochloride; Boehringer) allowed detection of cell nuclei. (Fig. 1a) Most of the cells expressed p53 histidine fusion protein clearly demonstrated by immunofluorescence with mAb 13/45/31 (1b) or rabbit anti-p53-serum (Fig 1c). With both antibodies the identical staining pattern of cell nuclei was observed. In another set of experiments we performed immunoblot analyses using mAb 13/45/31. Lysates of bacteria transformed with pQE/HBx comprising the X-reading frame of the hepatitis B virus and an enterokinase cleavage site cloned into pQE-8 were analysed before and after induction with IPTG (isopropyl- β -D-thiogalactopyranoside; Serva) by polyacrylamide gel electrophoresis (PAGE) and Coomassie-blue staining (Fig. 2, lanes 1 and 2).

Fig. 2 Detection of HBx protein with N-terminal histidine-tag and enterokinase cleavage site in bacterial lysates after PAGE and Coomassie-blue staining (lanes 1 and 2) and immunoblotting of a gel run in parallel with the mAb 13/45/31 (lanes 1' and 2'). Samples were analysed before (lanes 1 and 1') and after IPTG induction of the bacteria (lanes 2 and 2'). Marker, M; 10kDa ladder (Gibco-BRL).

Fig. 3 Detection of purified recombinant Cofilin expressed as fusion protein tagged at the N-terminus with an enterokinase site and six histidine residues. Immunoblotting was performed before (lanes 1 and 3) and after cleavage with enterokinase (lanes 2 and 4) with mAb 13/45/31 (lanes 1 and 2) or a polyclonal rabbit anti-Cofilin antibody (lanes 3 and 4). Marker positions are 15, 18, 30, 45 and 67 kDa, respectively.


Fig. 4 Detection of purified recombinant fusion proteins tagged at the C-terminus with six histidine residues by PAGE and Coomassie-blue staining (a) and immunoblotting (b) with mAb 13/45/31. Proteins were C22 (lanes 1 and 1'), C33 (lanes 2 and 2') and Ba (lanes 3 and 3'). Marker, M; 10 kDa ladder with

lowest band 20 kDa.

Fig. 5 Immunoprecipitation of recombinant histidine and enterokinase site tagged fusion protein hsp70 from bacterial extracts with mAb 13/45/31. Precipitation was performed with bacterial lysate before (lane 1) and after (lane 2) induction and analyzed following PAGE and immunoblotting with a polyclonal rabbit anti-hsp70 serum. Marker, M; 10kDa ladder.

The ~18kDa fusion protein was clearly detectable after IPTG induction (Fig. 2, lane 2). The corresponding immunoblot is shown in Fig. 2 (lane 1' and 2') and documents the high specificity of mAb 13/45/31. Even low amounts of protein not attributable to the recombinant protein after Coomassie-blue staining were detectable after immunoblotting with the mAb 13/45/31 (Fig. 2, lane 1 and 1'). Cleavage of the N-terminal-tagged histidine-enterokinase-fusion protein with enterokinase (Boehringer) prior to immunoblotting resulted in the complete loss of the specific signal of mAb 13/45/31 as shown for the Cofilin fusion protein (Fig. 3, compare lane 1 and 2). A blot prepared in parallel using a polyclonal rabbit anti-Cofilin serum served as a control (Fig. 3, lane 3 and 4). In a further assay we performed PAGE followed by immunoblotting with purified fusion proteins tagged solely with six histidine residues at the C-terminus. After Coomassie-blue staining (Fig. 4a) the pattern of the proteins C22 (lane 1 and 1'), C33 (lanes 2 and 2'), and Ba (lanes 3 and 3') was identical to that observed after immunoblotting using mAb 13/45/31 (Fig. 4b) thus demonstrating the highly specific affinity of mAb 13/45/31 to the epitop formed by the six histidine residues. To prove further applicabilities of mAb 13/45/31 we performed immunoprecipitation studies. One example is given in Fig. 5. Bacteria were transformed with a histidine-enterokinase-heat shock protein (hsp70) construct and expression was induced with IPTG. Using the mAb 13/45/31 we performed immunoprecipitation with bacterial lysates before and after induction. After PAGE of the precipitate and immunoblotting with a polyclonal rabbit anti-hp70 serum we detected the precipitated hsp70 fusion protein at the expected mol. wt of

~70 kDa. Other applications we have successfully examined are the use of mAb 13/45/31 in ELISA and microinjection studies (data not shown).

From the data presented here, using mAb 13/45/31 in immunofluorescence, immunoblotting and immunoprecipitation studies and, in addition, in analysis using ELISA and microinjection techniques, we believe mouse monoclonal anti-histidine tag antibody 13/45/31 is a powerful tool in molecular and cellular biology.

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