


Mouse Monoclonal Antibody (CAL2) Against All *CALRETICULIN (CALR)* Mutations

Product Information

Catalog No.:	DIA-CAL-250 (250µl) DIA-CAL-100 (100µl)	Practical implementation:	CAL2 labels the megakaryocytes in myeloproliferative neoplasms (essential thrombocythaemia (ET) and primary myelofibrosis (PMF)) with <i>CALR</i> mutation and enables to distinguish ET and PMF with <i>CALR</i> mutation from polycythemia vera (PV), from <i>CALR</i> mutation negative ET and PMF and from reactive bone marrow.
Clone:	CAL2	Staining pattern:	Cytoplasmic staining of megakaryocytes harboring <i>CALR</i> mutation. The CAL2 IHC assay indicates absence of <i>CALR</i> mutation when all megakaryocytes remain unlabeled. Formalin-fixed BM: Simple and fast analysis by a specific and intense staining of megakaryocytes and a variable amount of cells with small nucleus B5-fixed BM: Simple and fast analysis by a specific and intense staining of megakaryocytes. For analysis include only megakaryocytes. Sometimes weak non-specific staining of granulopoietic or erythropoietic cells
Isotype:	Mouse IgG2a		
Immunogen:	C-neotermius of mutated <i>CALR</i> .		
Specificity:	Human <i>CALR</i> protein expressed by all types of Exon 9 <i>CALR</i> mutations (deletion/insertion in 19p 13.3-13.2 of)		
Application:	Immunohistochemistry (IHC) for B5-fixed or formalin-fixed paraffin-embedded (FFPE) tissue with or without EDTA-decalcification		
Physical state:	Lyophilized powder		
Reagent provided:	Antibody purified from culture supernatant in PBS with 2% BSA, 0.05% NaN ₃ , pH 7.4.		
Storage and stability:	Unopened vials for at least 1 year at -20°C; As reconstituted liquid several weeks at 2-8 °C; For long term storage freeze at -20°C or -80°C. Avoid repeated freeze/thaw cycles.		
Instructions for use:	Reconstitute DIA-CAL-250 with 250 µl and DIA-CAL-100 with 100 µl sterile distilled water followed by gentle shaking for 10 min. Detection of <i>CALR</i> mut with clone CAL2, Heat-induced epitope retrieval (HIER): I. Formalin-fixed and EDTA-decalcified BM: Citrate-buffer pH6.0, 10min pressure cooker or 30-60min microwave at 98-100°C II. B5-fixed and EDTA-decalcified BM: EDTA-buffer pH8.0, 10min pressure cooker or 30-60min microwave at 98-100°C III. B5-fixed and acid-decalcified BM: EDTA-buffer pH8.0, 10min pressure cooker or EDTA-TRIS-buffer pH9.0, 60min microwave at 98-100°C After HIER follow standard IHC protocols. The CAL2 antibody IHC is suited for using automated platforms.		
Dilution:	1:20-1:40 for IHC	Positive control:	Megakaryocytes from <i>CALR</i> mutated PMNs
General recommendation:	Validation of antibody performance/protocol is the responsibility of the end user. Positive/negative controls should be run simultaneously with tissue specimen.	Negative control:	Megakaryocytes of reactive bone marrow specimens or <i>JAK2</i> mutated PV
		Safety notes:	The reconstituted liquid contains 0.05% sodium azide as a preservative. Avoid skin and eye contact, inhalation and ingestion.
		References:	1. Stein, H et al. A new monoclonal antibody (CAL2) detects <i>CALRETICULIN</i> mutations in formalin-fixed and paraffin embedded bone marrow sections. <i>Leukemia</i> 2015, accepted for publication. 2. Nangalia J et al. Somatic <i>CALR</i> Mutations in Myeloproliferative Neoplasms with Nonmutated <i>JAK2</i> . <i>N Engl J Med</i> 369(25): 2391-2405, 2013. 3. Klampfl T et al. Somatic Mutations of Calreticulin in Myeloproliferative Neoplasms. <i>N Engl J Med</i> 369(25): 2379-2390, 2013.



For research use only. Not for diagnostic or therapeutic use.

