



Instructions For Use LGA-IFU		
Rev. Date: July 20, 2008	Revision: 1	Page 1 of 1

P.O. Box 3286 - Logan, Utah 84323, U.S.A. - Tel. (800) 729-8350 - Fax (435) 755-0015 - www.scytek.com

Light Green Solution

Description: Light Green Solution is intended for use in histological applications as a general counterstain. When used correctly, various shades of green can be obtained to aid in visualization of tissue components.

Cytoplasm: Green
Nuclei: Green

Uses/Limitations: For In-Vitro Diagnostic use only.
Histological and Cytological applications.
Do not use past expiration date.
Use caution when handling these reagents.

Control Tissue: Any fixed or frozen tissue.
Cell smear.


Storage: Room Temperature (18-25° C.)


Availability:	<u>Item #</u>	<u>Volume</u>
	LGA125	125 ml
	LGA500	500 ml
	LGA999	1000 ml
	LGA-10000	10 Liters
	LGA-20000	20 Liters

Precautions: Avoid contact with skin and eyes.
Harmful if swallowed.
Follow all Federal, State, and local regulations regarding disposal.
Use in chemical fume hood whenever possible.
Wear protective clothing.

Procedure (Standard):

1. Deparaffinize sections if necessary and hydrate to distilled water.
2. Rinse slide in distilled water.
3. Stain slide in Light Green Solution for 2 minutes.
4. Dehydrate through graded alcohols.
5. Clear, and mount in synthetic resin.

Storage: 18° C  25° C

 ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
435-755-9848
U.S.A.

 Authorized Representative in Europe
 (Regulatory affairs only)
EmergoEurope (31)(0) 70 345-8570
Molsnstraat 15
2513 BH Hague, The Netherlands