



## Instructions For Use LGD-IFU

Rev. Date: Oct. 18, 2011

Revision: 1

Page 1 of 1

P.O. Box 3286 - Logan, Utah 84323, U.S.A. - Tel. (800) 729-8350 – Tel. (435) 755-9848 - Fax (435) 755-0015 - [www.scytek.com](http://www.scytek.com)

# Light Green Solution (0.1%)

### Description:

Light Green Solution (0.1%) is intended for use in histological applications as a counterstain for Grocott's Silver Stain in addition to other procedures. When used correctly, various shades of green can be obtained to aid in visualization of tissue components.

Cytoplasm: Light Green  
Nuclei: Light Green/Blue

### Uses/Limitations:

Not to be taken internally.  
For In-Vitro Diagnostic use only.  
Histological applications.  
Do not use past expiration date.  
Use caution when handling reagents.  
Non-Sterile.

### Availability:

<u>Item #</u>	<u>Volume</u>
LGD125	125 ml
LGD500	500 ml
LGD999	1000 ml

### Storage:

Store at 18-25°C.

### Precautions:

Avoid contact with skin and eyes.  
Harmful if swallowed.  
Follow all Federal, State, and local regulations regarding disposal.

### Procedure:

1. Deparaffinize sections if necessary and hydrate to distilled water.
2. Rinse slide in distilled water.
3. Stain slide in Light Green Solution (0.1%) for 1-2 minutes.
4. Rinse 1 time in Distilled/Deionized water.
5. Dehydrate quickly through graded alcohols.
6. Clear, and mount in synthetic resin.

Storage: 18° C



25° C



ScyTek Laboratories, Inc.  
205 South 600 West  
Logan, UT 84321  
U.S.A.



 EmergoEurope (31)(0) 70 345-8570  
Molsnstraat 15  
2513 BH Hague, The Netherlands