



Instructions For Use ORG-IFU		
Rev. Date: July 14, 2009	Revision: 1	Page 1 of 2

P.O. Box 3286 - Logan, Utah 84323, U.S.A. - Tel. (800) 729-8350 - Fax (435) 755-0015 - www.scytek.com

Oil Red O Solution (For Fat)

Description: Oil Red O Solution is a component of the Oil Red O Stain Kit (For Fat) (Catalog# ORK-1) intended for use in the histological visualization of fat cells and neutral fat. This kit may be used **ONLY** on frozen tissue sections, fresh smears, or touch preps.

Fat Cells:	Red
Neutral Fat:	Red
Nuclei:	Blue

Uses/Limitations: For In-Vitro Diagnostic use only. Histological applications. Do not use past expiration date.

Control Tissue: Any frozen section containing fat.

Availability/Contents:

<u>Item #</u>	<u>Kit Contents</u>	<u>Volume</u>	<u>Storage</u>
ORG125	Oil Red O Solution	125 ml	Room Temperature
ORG500	Oil Red O Solution	500 ml	Room Temperature
ORG999	Oil Red O Solution	1000 ml	Room Temperature

Required But Not Included:


<u>Item #</u>	<u>Kit Contents</u>	<u>Volume</u>	<u>Storage</u>
PRG500	Propylene Glycol	500 ml	Room Temperature
HMM125	Hematoxylin, Mayer's (Lillie's Mod.)	125 ml	Room Temperature

Precautions: Keep away from open flame. Avoid contact with skin and eyes. Harmful if swallowed. Follow all Federal, State, and local regulations regarding disposal. Use in chemical fume hood whenever possible.

Procedure (Standard):

1. Prepare fresh or frozen tissue section as usual.
 2. Place slide in Propylene Glycol for 2 minutes.
 3. Incubate slide in Oil Red O Solution for 6 minutes.
- Note: Prepare mixture of 85% Propylene Glycol in distilled water.
4. Differentiate tissue section in 85% Propylene Glycol for 1 minute.

Storage: 18° C  25° C

 ScyTek Laboratories, Inc.
 205 South 600 West
 Logan, UT 84321
 435-755-9848
 U.S.A.

 
 (Regulatory affairs only)
 EmergoEurope (31)(0) 70 345-8570
 Molsstraat 15
 2513 BH Hague, The Netherlands

Instructions For Use ORG-IFU

Rev. Date: July 14, 2009

Revision: 1

Page 2 of 2

P.O. Box 3286 - Logan, Utah 84323, U.S.A. - Tel. (800) 729-8350 - Fax (435) 755-0015 - www.scytek.com

5. Rinse slide in 2 changes of distilled water.
6. Incubate slide in Hematoxylin, Mayer's (Lillie's Modification) for 1-2 minutes.
7. Rinse slide thoroughly in tap water.
8. Rinse slide in 2 changes of distilled water.
9. Coverslip using an aqueous mounting medium (cat# AML060).

References:

1. Sheehan, DC., Hrapchak, BB. Theory and Practice of Histotechnology; 1980, page 225.

Storage: 18° C




25° C



ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
435-755-9848
U.S.A.



 (Regulatory affairs only)
EmergoEurope (31)(0) 70 345-8570
Molsstraat 15
2513 BH Hague, The Netherlands