



DATA SHEET

CD79a/mb-1 (B Cell Marker) Ab-1 (Clone HM47/A9)

Mouse Monoclonal Antibody

Cat. #DLN-09117, DLN-09118, or DLN-09116 (0.1ml, 0.5ml, or 1.0ml at 200µg/ml) (Purified Ab with BSA and Azide)

Cat. #DLN-09119 or DLN-09120 (0.1ml or 0.2ml at 1.0mg/ml) (Purified Ab without BSA and Azide)

Cat. #DLN-09114, DLN-09115, or DLN-09113 (0.1ml, 0.5ml, or 1.0ml at 200µg/ml) (Biotin-Labeled Ab with BSA and Azide)

Cat. #DLN-09121 (7.0ml) (Ready-to-Use for Immunohistochemical Staining)

Description: A disulphide-linked heterodimer, consisting of mb-1 (or CD79a) and B29 (or CD79b) polypeptides, is non-covalently associated with membrane-bound immunoglobulins on B cells to constitute the B cell Ag receptor. CD79a first appears at pre B cell stage and persists until the plasma cell stage where it is found as an intracellular component. CD79a is found in the majority of acute leukemias of precursor B cell type, in B cell lines, B cell lymphomas, and in some myelomas. It is not present in myeloid or T cell lines.

Mol. Wt. of Antigen: 47kDa

Epitope: GTYQDVGSLNIADVQ

Species Reactivity: Human¹, Monkey¹, Horse¹, Cow¹, Pig¹, Rabbit¹, Guinea pig¹, Rat¹, Mouse¹ and Chicken (Plasma cells). ¹ Does not react with opossum. ¹

Clone Designation: HM47/A9

Ig Isotype / Light Chain: IgG₁ / κ

Immunogen: A synthetic peptide consisting of GTYQDVGSLNIADVQ of human CD79a protein.

Applications and Suggested Dilutions:

- Immunofluorescence
 - Western Blotting (0.5-1µg/ml for 2hrs at RT))
 - Immunohistology (Formalin/paraffin)
(Ab 0.5-1.0µg/ml for 30 min at RT)
- * [Staining of formalin-fixed tissues REQUIRES boiling tissue sections in 10mM citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.]

The optimal dilution for a specific application should be determined by the investigator.

Positive Control: Raji cells. Lymph node or tonsil.

Cellular Localization: Cell membrane

Supplied As: 200µg/ml Ab purified from ascites fluid by Protein G chromatography. Prepared in 10mM PBS, pH 7.4, with 0.2% BSA and 0.09% sodium azide. Also available without BSA and azide at 1mg/ml, **or** Prediluted antibody which is ready-to-use for staining of formalin-fixed, paraffin-embedded tissues.

dianova GmbH

Warburgstr. 45 • 20354 Hamburg

Telefon (040)45067-0 • Telefax (040) 45067-490 • www.dianova.de



DATA SHEET

CD79a/mb-1 (B Cell Marker) Ab-1 (Clone HM47/A9)

Mouse Monoclonal Antibody

Cat. #DLN-09117, DLN-09118, or DLN-09116 (0.1ml, 0.5ml, or 1.0ml at 200µg/ml) (Purified Ab with BSA and Azide)

Cat. #DLN-09119 or DLN-09120 (0.1ml or 0.2ml at 1.0mg/ml) (Purified Ab without BSA and Azide)

Cat. #DLN-09114, DLN-09115, or DLN-09113 (0.1ml, 0.5ml, or 1.0ml at 200µg/ml) (Biotin-Labeled Ab with BSA and Azide)

Cat. #DLN-09121 (7.0ml) (Ready-to-Use for Immunohistochemical Staining)

Storage and Stability:

Ab with sodium azide is stable for 24 months when stored at 2-8°C. Antibody WITHOUT sodium azide is stable for 36 months when stored at below 0°C.

Key References:

1. Jones M, et. al. (1993) J Immunol, 150:5429-5435.
2. Mason D Y, et. al. (1991) J Immunol 147:2474-2482.

Limitations and Warranty:

Our products are intended FOR RESEARCH USE ONLY and are not approved for clinical diagnosis, drug use or therapeutic procedures. No products are to be construed as a recommendation for use in violation of any patents. We make no representations, warranties or assurances as to the accuracy or completeness of information provided on our data sheets and website. Our warranty is limited to the actual price paid for the product. Dianova is not liable for any property damage, personal injury, time or effort or economic loss caused by our products.

Material Safety Data:

This product is not licensed or approved for administration to humans or to animals other than the experimental animals. Standard Laboratory Practices should be followed when handling this material. The chemical, physical, and toxicological properties of this material have not been thoroughly investigated. Appropriate measures should be taken to avoid skin and eye contact, inhalation, and ingestion. The material contains 0.09% sodium azide as a preservative. Although the quantity of azide is very small, appropriate care should be taken when handling this material as indicated above. The National Institute of Occupational Safety and Health has issued a bulletin citing the potential explosion hazard due to the reaction of sodium azide with copper, lead, brass, or solder in the plumbing systems. Sodium azide forms hydrazoic acid in acidic conditions and should be discarded in a large volume of running water to avoid deposits forming in metal drainage pipes.

For Research Use Only

dianova GmbH

Warburgstr. 45 • 20354 Hamburg

Telefon (040)45067-0 • Telefax (040) 45067-490 • www.dianova.de