

## Product Information

|                                   |   |
|-----------------------------------|---|
| <b>Catalog Number:</b>            | 009-000-050   |
| <b>Product:</b>                   | ChromPure Human Transferrin   |
| <b>Physical State:</b>            | Sterile-filtered liquid   |
| <b>Size:</b>                      | 5 mg  |
| <b>Transferrin Concentration:</b> | 5.9 mg/ml   |
| <b>Conjugation:</b>               | unconjugated  |
| <b>Buffer:</b>                    | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6  |
| <b>Preservative:</b>              | None  |
| <b>Purity:</b>                    | Based on immunoelectrophoresis at an antigen concentration of 20 mg/ml, the pattern of precipitation against anti-human whole serum is the same as that against anti-human transferrin. No precipitin line was detected against anti-human IgG (H+L). |
| <b>Storage:</b>                   | Store at 2-8°C under sterile conditions. Prepare working dilution fresh each day. Expiration date: one year from date of receipt. The expiration date may be extended if test results are acceptable for the intended use.                            |
| <b>Safety Information:</b>        | This product contains no substances considered to be hazardous to health. Although the product is classified as non-hazardous, Standard Laboratory Practices should be followed when handling this material.  |

**Note:** This information sheet contains general product information only. For detailed lot specific information consult the vial label and the data sheet supplied by the manufacturer upon delivery of the product.

**This product is for In Vitro experimental use only. Not for Therapeutic or Diagnostic use.**

### Manufacturer:

**Jackson ImmunoResearch Laboratories, Inc.**  
872 West Baltimore Pike  
West Grove, PA, USA 19390

### Distribution and Support:

**dianova GmbH**  
Warburgstraße 45  
20354 Hamburg  
www.dianova.com

**Technical Support:**  
Phone: 040 45 06 70  
Email: [info@dianova.de](mailto:info@dianova.de)

**Ordering Information:**  
Phone: 040 45 06 70  
Fax: 040 45 06 74 90  
Email: [order@dianova.de](mailto:order@dianova.de)