



Product Information

Catalog Number:	109-125-011
Product:	PerCP-AffiniPure Goat Anti-Human IgA, α Chain Specific
Physical State:	Freeze-dried solid
Size:	0.5 ml
Antibody Concentration:	± 0.5 mg/ml (exact concentration lot dependent)
Fluorophore:	Peridinin-Chlorophyll Protein Complex, purified from Dinophyceae sp. Amax = 482 nm; Emax = 677 nm
Suggested Dilution Range:	1:25 - 1:100 for most applications
Buffer:	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6
Stabilizer:	15 mg/ml Bovine Serum Albumin (IgG-Free, Protease-Free)
Preservative:	0.05% Sodium Azide
Purity:	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Antibody Specificity:	Based on immunoelectrophoresis and/or ELISA, the antibody reacts with the heavy chain of human IgA but not with human IgG or IgM. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with IgA from other species.
Reconstitution and Storage:	Store freeze-dried solid at 2-8°C. Rehydrate with 0.5 ml dH2O and centrifuge if not clear. Store at 2-8°C - do not freeze. Prepare working dilution on day of use. Expiration date:six months from date of rehydration. The expiration date may be extended if test results are acceptable for the intended use.
Safety Information:	Standard Laboratory Practices should be followed when handling this material. The chemical, physical, and toxicological properties of this material have not been thoroughly investigated. The antibody contains 0.05% sodium azide (a poisonous and hazardous substance) as preservative. Although this concentration is not regarded as dangerous to health, appropriate measures should be taken to avoid skin and eye contact, inhalation, and ingestion.

Note: This information sheet contains general product information only. For detailed lot specific information consult the vial label and the data sheet supplied by the manufacturer upon delivery of the product.

This product is for In Vitro experimental use only. Not for Therapeutic or Diagnostic use.

Manufacturer:

Jackson ImmunoResearch Laboratories, Inc. 872 West Baltimore Pike West Grove, PA, USA 19390

Distribution and Support:

dianova GmbH Warburgstraße 45 20354 Hamburg www.dianova.com Technical Support: Phone: 040 45 06 70 Email: info@dianova.de **Ordering Information:** Phone: 040 45 06 70 Fax: 040 45 06 74 90 Email: <u>order@dianova.de</u>