



## **Product Information**

Catalog Number:	200-622-156
Product:	Alexa Fluor® 680-conjugated IgG Fraction Monoclonal Mouse Anti-Digoxin
Physical State:	Freeze-dried solid
Size:	0.3 mg (This product has been overfilled to ensure total recovery of stated quantity.)
IgG Concentration:	1.7 mg/ml
Fluorophore:	Alexa Fluor®680 carboxylic acid Amax = 684 nm; Emax = 700 nm
Suggested Dilution Range:	1:50,000 - 1:200,000 for Western blotting
Buffer:	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6
Stabilizer:	15 mg/ml Bovine Serum Albumin (IgG-Free, Protease-Free)
Preservative:	0.05% Sodium Azide
Purity:	The IgG fraction was purified from mouse ascites by chromatographic procedures not involving
	affinity purification.
Reconstitution and Storage:	affinity purification. Store freeze-dried powder at 2-8°C. When ready to use, rehydrate with 0.2 ml dH2O and centrifuge if not clear. Product is stable for about 6 weeks at 2-8°C as an undiluted liquid. Prepare working dilution fresh each day. For extended storage after rehydration, aliquot and freeze at -70°C or below. Avoid repeated freezing and thawing. Alternatively, add an equal volume of glycerol (ACS grade or better) for a final concentration of 50%, and store at -20°C as a liquid. Note: adding glycerol reduces the stated protein concentration and dilution range by one-half. Expiration date: one year from date of rehydration. The expiration date may be extended if test results are acceptable for the intended use.

**Note:** This information sheet contains general product information only. For detailed lot specific information consult the vial label and the data sheet supplied by the manufacturer upon delivery of the product.

This product is for In Vitro experimental use only. Not for Therapeutic or Diagnostic use.

Manufacturer:

Jackson ImmunoResearch Laboratories, Inc. 872 West Baltimore Pike West Grove, PA, USA 19390

**Distribution and Support:** 

dianova GmbH Warburgstraße 45 20354 Hamburg www.dianova.com Technical Support: Phone: 040 45 06 70 Email: info@dianova.de Ordering Information: Phone: 040 45 06 70 Fax: 040 45 06 74 90 Email: order@dianova.de