



Rabbit F(ab')2 anti-Chicken IgG (H+L)-unconj., MinX none

General Information

Catalog Number: SEC-182655	Physical State:	Lyophilized			
Quantity: 20 mg	Reconstitution Volume:	2.0 mL			
	Reconstitution Buffer:	Restore with deionized water (or equivalent) at room temperature			
	Shipping Conditions:				
	Product Expiration	Expiration date is one (1) year from date of opening.			
Antibody Host / Format	I				
Host Species: Rabbit	Clonality / Isotype / Clor	ne: polyclonal Ig			
Antibody Format: IgG F(ab')2 Fragment	Concentration:	12.5 mg/mL by UV absorbance at 280 nm			
	Conjugation:	unconjugated			
	Maximum Excitation / E	mission: /			
Specificity	 				
Target Species: Chicken	Immunogen:	Chicken IgG whole molecule			
Antibody Specificity: IgG (H+L)	This product is a F(ab')2 fragment of IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation, ion exchange chromatography and pepsin digestion followed by chromatographic separation and extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Rabbit Serum, Chicken IgG and Chicken Serum. No reaction was observed against anti-Rabbit IgG F(c) or anti-Pepsin.				
	Minimal Cross Reactivity: none				
Formulation, Transport and Storage	I				
Storage Buffer:0.01 M Sodium Phosphate,0.15 M Sodium Chloride, pH 7.2Stabilizer:NonePreservatives:0.01% (w/v) Sodium Azide	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.				
	Storage Temperature: 4-8°C				
Application Recommendation	I				

This product is for In Vitro exp	erimental use only. Not for Ther	apeutic or Diagnostic use.				
Manufacturer				Manag	gement m	
BIOZOL Diagnostica Vertrieb GmbH Leipziger Straße 4	Phone +49 (89) 3799 666 6 Fax +49 (89) 3799 666 99 E-mail info@biozol.de	www.dianova.com Your Secondary Antibody Portal	TÜVRhein	R	001:2015	
85386 Eching	www.biozol.de	Powered by BIOZOL	ZERTIFIZI	ERT www.tu ID 9000		



Product Datasheet / Instruction for use



- **Application:** ELISA, Immunohistochemistry (IHC), Immunohistochemistry (Paraffin-embedded Sections), Immunohistochemistry (frozen sections), Western Blot
- Dilution: ELISA 1:20.000 1:100.000, Immunohistochemistry (IHC) 1:1,000 1:5,000, Western Blot (WB) 1:2.000 1:10.000

Application Note:

Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity. The maximum amount of reagent required to stain 1 x 10E6 cells in flow cytometry is approximately 1.0 µg of antibody. Lesser amounts of reagent may be sufficient for staining. Optimal titers for other applications should be determined by the researcher. As a general guideline dilutions of 1:100 to 1:250 should be suitable for most applications.

Background Information

F(ab')2 Antibody was generated by enzymatic cleavage and subsequent separation from the Fc fragment. Because of their smaller size, F(ab)2 fragments offer several advantages over intact antibodies for use in certain immunochemical techniques and experimental applications. F(ab)2 fragments penetrate into tissue samples and show better antigen recognition and signal generation in IHC. F(ab)2 fragments lack the Fc region and therefore do not bind Fc receptors which effectively lowers background staining. F(ab')2 Antibody is ideal for investigators who routinely perform flow cytometry, immunohistochemistry or IHC and other immunoassays.

Safety Information

This reagent contains sodium azide (a poisonous and hazardous substance) as preservative. Although this concentration is not regarded as dangerous to health, appropriate measures should be taken to avoid skin and eye contact, inhalation, and ingestion. Standard Laboratory Practices should be followed when handling this material.

Disclaimer

This product is warranted to perform in conformance with product specifications and to be free from defects in material and workmanship. Products are supplied for research use only. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated, this warranty is limited to one year from date of shipment when the product is subjected to normal, proper and intended usage. Buyer's exclusive remedy for non-conforming products during the warranty period is limited to replacement of or refund for the non-conforming product(s). There is no obligation to replace products as the result of accident or disaster or inappropriate use of the products in a manner for which they were not designed, or improper storage and handling of the products.

This product is for In Vitro experimental use only. Not for Therapeutic or Diagnostic use.

Manufacturer

BIOZOL Diagnostica Vertrieb GmbH Leipziger Straße 4 85386 Eching Phone +49 (89) 3799 666 6 Fax +49 (89) 3799 666 99 E-mail info@biozol.de www.biozol.de

www.dianova.com Your Secondary Antibody Portal

Powered by BIOZOL

TÜVRheinland ZERTIFIZIERT

Management System ISO 9001:2015