

Product Datasheet / Instruction for use



# Rabbit F(ab')2 anti-Horse IgG (H+L)-HRPO, MinX none

## **General Information**

Catalog Number:	<b>bg Number:</b> SEC-183714 <b>Physical State:</b> Lyophilized				
Quantity:	0,5 mg	Reconstitution Volume:	500 μL		
		<b>Reconstitution Buffer:</b>	Restore with deionized water (or equivalent)		
		Shipping Conditions:	at room temperature		
		Product Expiration Expiration date is one (1) opening.			
Antibody Host / Fo	ormat				
Host Species:	Rabbit	Clonality / Isotype / Clon	ne: polyclonal Ig		
Antibody Format:	IgG F(ab')2 Fragment	Concentration:	1.0 mg/mL by UV absorbance at 280 nr		
		Conjugation:	Horseradish Peroxidase (HRP)		
		Maximum Excitation / Emission: /			
Specificity					
Target Species:	Horse	Immunogen:	Horse IgG whole molecule		
Antibody Specificity: IgG (H+L)		This product was prepared from monospecific antiserum by immunoaffinity chromatography using Horse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation. Assay b immunoelectrophoresis resulted in a single precipitin arc against anti- Peroxidase, anti-Rabbit Serum, Horse IgG and Horse Serum. No reactio was observed against anti-Pepsin or anti-Rabbit IgG F(c).			
		Minimal Cross Reactivity: none			
Formulation, Trans	sport and Storage				
Storage Buffer:0.02 M PotassiumPhosphate,0.15 M Sodium Chloride,Phosphate,0.15 M Sodium Chloride,Phosphate,10 mg/mL Bovine SerumAlbumin (BSA) - Immunoglobulin and Proteasefree		Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.			
Preservatives: 0.01% (w/v) Gentamicin Sulfate		Storage Temperature: 4-8°C			
Application Recom	imendation				

This product is for In Vitro experimental use only. Not for Therapeutic or Diagnostic use.							
Manufacturer				Management System			
BIOZOL Diagnostica Vertrieb GmbH Leipziger Straße 4	Phone +49 (89) 3799 666 6 Fax +49 (89) 3799 666 99 E-mail info@biozol.de	www.dianova.com	TÜVRheinland	ISO 9001:2015			
85386 Eching	www.biozol.de	Powered by BIOZOL	ZERTIFIZIERT	www.tuv.com ID 9000019771			



#### Product Datasheet / Instruction for use



Application: ELISA, Immunohistochemistry (IHC), Immunohistochemistry (Paraffin-embedded Sections), Immunohistochemistry (frozen sections), Western Blot

**Dilution:** ELISA 1:90.000, Immunohistochemistry (IHC) 1:500 - 1:2,500, Western Blot (WB) 1:1.000 - 1:5.000

#### **Application Note:**

This product is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.

### **Background Information**

F(ab')2 Anti-Horse IgG Peroxidase Antibody was generated by enzymatic cleavage and subsequent separation from the Fc fragment. Because of their smaller size, F(ab)2 fragments offer several advantages over intact antibodies for use in certain immunochemical techniques and experimental applications. F(ab)2 fragments penetrate tissue samples and show better antigen recognition and signal generation in IHC. F(ab)2 fragments lack the Fc region and therefore do not bind Fc receptors which effectively lowers background staining. F(ab')2 Antibody is ideal for investigators who routinely perform flow cytometry, immunohistochemistry or IHC and other immunoassays.

### **Safety Information**

This reagent contains Gentamicin Sulfate (a poisonous and hazardous substance) as preservative. Although this concentration is not regarded as dangerous to health, appropriate measures should be taken to avoid skin and eye contact, inhalation, and ingestion. Standard Laboratory Practices should be followed when handling this material.DO NOT ADD SODIUM AZIDE! Use of Sodium Azide as a preservative will substantially inhibit the enzyme activity of horseradish peroxidase!

#### Disclaimer

This product is warranted to perform in conformance with product specifications and to be free from defects in material and workmanship. Products are supplied for research use only. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated, this warranty is limited to one year from date of shipment when the product is subjected to normal, proper and intended usage. Buyer's exclusive remedy for non-conforming products during the warranty period is limited to replacement of or refund for the non-conforming product(s). There is no obligation to replace products as the result of accident or disaster or inappropriate use of the products in a manner for which they were not designed, or improper storage and handling of the products.

This product is for In Vitro experimental use only. Not for Therapeutic or Diagnostic use.

Manufacturer

BIOZOL **Diagnostica Vertrieb GmbH** Leipziger Straße 4 85386 Eching

Phone +49 (89) 3799 666 6 Fax +49 (89) 3799 666 99 E-mail info@biozol.de www.biozol.de

www.dianova.com Your Secondary Antibody Portal

Powered by BIOZOL



Management System ISO 9001:2015