

BIOZOL FIT FOR SCIENCE SEC-183732 2023-02-23 1/2

Goat F(ab')2 anti-Human Kappa light chain-FITC, MinX none

General Information

Catalog Number: SEC-183732	Physical State:	Lyophilized		
Quantity: 1 mg	Reconstitution Volume:	1.0 mL		
	Reconstitution Buffer:	Restore with deionized water (or equivalent) at room temperature Expiration date is one (1) year from date of opening.		
	Shipping Conditions:			
	Product Expiration			
Antibody Host / Format	l			
Host Species: Goat	Clonality / Isotype / Clone: polyclonal Ig			
Antibody Format: IgG F(ab')2 Fragment	Concentration:	1.0 mg/ml by UV absorbance at 280 nm		
	Conjugation:	FITC (Fluoresceinisothiocyanat)		
	Maximum Excitation / Emission: 492 nm / 520 nm			
Specificity	I			
Target Species: Human	Immunogen:	Human kappa light chain		
Antibody Specificity: Kappa light chain	F(ab')2 HUMAN kappa Antibody Fluorescein Conjugated was prepared from monospecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-fluorescein, anti-Goat Serum, Huma IgG and Human Serum. No reaction was observed against anti-Pepsin or anti-Goat IgG F(c). Specificity was confirmed by ELISA at less than 1% cross reactivity against other human heavy or light chain isotypes.			
	Minimal Cross Reactivity: none			
Formulation, Transport and Storage				
Storage Buffer:0.02 M PotassiumPhosphate, 0.15 M Sodium Chloride, pH 7.2Stabilizer:10 mg/mL Bovine SerumAlbumin (BSA) - Immunoglobulin and Protease	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.			
free				
Preservatives: 0.01% (w/v) Sodium Azide	Storage Temperature:	4-8°C		

This product is for In Vitro experimental use only. Not for Therapeutic or Diagnostic use.					
Manufacturer				Management System	
BIOZOL Diagnostica Vertrieb GmbH Leipziger Straße 4	Phone +49 (89) 3799 666 6 Fax +49 (89) 3799 666 99 E-mail info@biozol.de	www.dianova.com Your Secondary Antibody Portal	TÜVRheinland	ISO 9001:2015	
85386 Eching	www.biozol.de	Powered by BIOZOL	ZERTIFIZIERT	www.tuv.com ID 9000019771	



Product Datasheet / Instruction for use



Application: FLISA, Flow Cytometry, Immunofluorescence

Dilution: FLISA 1:10.000 - 1:50.000, Flow cytometry 1:500 - 1:2,500, Immunofluorescence 1:1,000 - 1:5,000

Application Note:

Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.

Background Information

F(ab')2 HUMAN kappa Antibody Fluorescein Conjugated is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.

Safety Information

This reagent contains sodium azide (a poisonous and hazardous substance) as preservative. Although this concentration is not regarded as dangerous to health, appropriate measures should be taken to avoid skin and eye contact, inhalation, and ingestion. Standard Laboratory Practices should be followed when handling this material.

Disclaimer

This product is warranted to perform in conformance with product specifications and to be free from defects in material and workmanship. Products are supplied for research use only. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated, this warranty is limited to one year from date of shipment when the product is subjected to normal, proper and intended usage. Buyer's exclusive remedy for non-conforming products during the warranty period is limited to replacement of or refund for the non-conforming product(s). There is no obligation to replace products as the result of accident or disaster or inappropriate use of the products in a manner for which they were not designed, or improper storage and handling of the products.

This product is for In Vitro experimental use only. Not for Therapeutic or Diagnostic use.

Manufacturer

BIOZOL **Diagnostica Vertrieb GmbH** Leipziger Straße 4 85386 Eching

Phone +49 (89) 3799 666 6 Fax +49 (89) 3799 666 99 E-mail info@biozol.de www.biozol.de

www.dianova.com Your Secondary Antibody Portal

Powered by BIOZOL



Management System ISO 9001:2015