

ScyTek
Laboratories

Instructions For Use
FGN-IFU

Rev. Date: 10/28/03

Revision: 2

Page 1 of 1

P.O. Box 3286 - Logan, Utah 84323, U.S.A. - Tel. (800) 729-8350 - Fax (435) 755-0015 - www.scytek.com

Fast Green FCF (Aqueous)

Description: Fast Green FCF is a stain with histological applications.

Availability:

<u>Item #</u>	<u>Volume</u>
FGN500	500 ml
FGN999	1000 ml
FGN010	10 Liters

Procedure:

1. Re-hydrate tissue as usual.
2. Stain cells (muscle) 5 minutes with Fast Green FCF.
3. Differentiate in 1% acetic acid. (AAE999)
4. Rinse in 95% alcohol 5 min.
5. Dehydrate, clear, and mount in permanent mounting medium.