

ScyTek
Laboratories

Instructions For Use SOH-IFU		
Rev. Date: May 4, 2007	Revision: 1	Page 1 of 1

P.O. Box 3286 - Logan, Utah 84323, U.S.A. - Tel. (800) 729-8350 - Fax (435) 755-0015 - www.scytek.com

Safranin O Solution


Description: Safranin O is commonly used for counterstaining nuclei red.


Availability:

<u>Item #</u>	<u>Volume</u>
SOH500	500 ml
SOH999	1000 ml

Procedure:

1. Re-hydrate tissue as usual.
2. Stain tissue section for 5 minutes with Safranin O Solution.
3. Rinse in 95% alcohol 1 minute followed by 1 minute in absolute alcohol.
4. Clear in xylene, and mount in permanent mounting medium.

Storage: 18° C  25° C

 ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
435-755-9848
U.S.A.

 Authorized Representative in Europe
 (Regulatory affairs only)
EmergoEurope (31)(0) 70 345-8570
Molsnstraat 15
2513 BH Hague, The Netherlands