



Goat Antiserum anti-Bovine IgG (H+L)-unconj., MinX none

General Information

Catalog Number: SEC-182407	Physical State:	Lyophilized			
Quantity: 2 ml	Reconstitution Volume:	2.0 mL			
	Reconstitution Buffer:	Restore with deionized water (or equivalent)			
	Shipping Conditions:	at room temperature			
	Product Expiration	Expiration date is one (1) year from date of opening.			
Antibody Host / Format					
Host Species: Goat	Clonality / Isotype / Clor	e: polyclonal Ig			
Antibody Format: Antiserum	Concentration:	28.5 mg/mL by Refractometry			
	Conjugation:	unconjugated			
	Maximum Excitation / E	mission: /			
Specificity					
Target Species: Bovine	Immunogen:	Bovine IgG whole molecule			
Antibody Specificity: IgG (H+L)	This product was prepared from monospecific antiserum by a delipidati and defibrination. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-goat serum, Bovine IgG and Bovine Serum.				
	Minimal Cross Reactivity: none				
Formulation, Transport and Storage	1				
Storage Buffer: 0.01 M Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and				
Stabilizer: None	thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an				
Preservatives: None	undiluted liquid. Dilute c	only prior to immediate use.			
	Storage Temperature:	4-8°C			
Application Recommendation					

Application: ELISA, Immunohistochemistry (IHC), Immunohistochemistry (Paraffin-embedded Sections), Immunohistochemistry (frozen sections)

Dilution: ELISA 1:20.000 - 1:100.000, Immunohistochemistry (IHC) 1:1,000 - 1:5,000, Western Blot (WB) 1:2.000 - 1:10.000

This product is for In Vitro experimental use only. Not for Therapeutic or Diagnostic use.						
Manufacturer				Management System		
BIOZOL Diagnostica Vertrieb GmbH Leipziger Straße 4	Phone +49 (89) 3799 666 6 Fax +49 (89) 3799 666 99 E-mail info@biozol.de	www.dianova.com Your Secondary Antibody Portal	TÜVRheinland	ISO 9001:2015		
85386 Eching	www.biozol.de	Powered by BIOZOL	ZERTIFIZIERT	www.tuv.com ID 9000019771		



Product Datasheet / Instruction for use



Application Note:

Anti-Bovine IgG Secondary Antibodies are designed for Western Blotting, ELISA and Immunohistochemistry. Secondary antibodies can also be used for a variety of other applications such as Assay Development, Flow Cytometry and in situ hybridization.

Background Information

Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G binds to viruses, bacteria, as well as fungi and facilitates their destruction or neutralization via agglutination (and thereby immobilizing them), activation of the compliment cascade, and opsinization for phagocytosis. This product possesses the F(ab) region possessing the epitope-recognition site, both heavy and light chains of the antibody molecule are present.

Safety Information

The product contains no hazardous constituents. Standard Laboratory Practices should be followed when handling this material.

Disclaimer

This product is warranted to perform in conformance with product specifications and to be free from defects in material and workmanship. Products are supplied for research use only. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated, this warranty is limited to one year from date of shipment when the product is subjected to normal, proper and intended usage. Buyer's exclusive remedy for non-conforming products during the warranty period is limited to replacement of or refund for the non-conforming product(s). There is no obligation to replace products as the result of accident or disaster or inappropriate use of the products in a manner for which they were not designed, or improper storage and handling of the products.

This product is for In Vitro experimental use only. Not for Therapeutic or Diagnostic use.

Manufacturer BIOZOL

Diagnostica Vertrieb GmbH Leipziger Straße 4 85386 Eching

Phone +49 (89) 3799 666 6 Fax +49 (89) 3799 666 99 E-mail info@biozol.de www.biozol.de

www.dianova.com Your Secondary Antibody Portal

Powered by BIOZOL



Management System ISO 9001:2015